# Why are these kits important?

In the event of a traumatic injury, these kits give a person most of the tools they need to control the bleeding until help arrives. This can be the difference between life and death, or the loss of a limb.

A person who is bleeding severely can die from blood loss in just minutes.

# A Bleed Control Be Ready!

kit is an important item to train and equip private citizens to use these medically proven bleeding control practices to help save the lives of people who have suffered a traumatic injury.

# A Bleed control Be Ready!

kit is specifically designed to stop heavy bleeding and handle lifethreatening injuries caused by gunshot, stabbing, crushing and other traumatic injuries.



#### Wisconsin based company

Guardian CPR, LLC PO Box 124 Wautoma, WI 54982 920-229-4047

#### Working with:

Police

Fire

**EMS** 

Churches

**School Districts** 

Government Agencies
Private businesses

#### **Services Offered:**

AED- Sales/Service
CPR Training
First Aid Training
Customizable First Aid Kits

# Bleed Control E Ready &









guardiancprllc.com 920-229-4047 With our **Bleed Control Be**Ready! kits, they will help
decrease the number of
preventable deaths due to
hemorrhaging with proper training
and use

We have four standard kits which include the following. Kits include:

Kit # 1 Contents:

- 1- Tourniquet
- 1- Permanent Marker
- 1- Trauma Shears
- 2- Pair Large Nitrile Gloves
- 1- Combat Application Tourniquet
- 2-5 x 9 Sterile Dressing
- 2-4" x 12.3' Conforming Stretch Gauze
- 1- Large Mylar Blanket

Kit # 2 add all from above plus 1- Israeli Type Bandage

Kit # 3 add all from above plus 2- Chest seals

Kit # 4 add all from above plus 1- Z-Fold Compress Gauze



### Kits can be customized

Need something not listed?

Just Ask!

# How best to control severe bleeding?

5 methods to control severe bleeding

- \* Call 911
- \* Place a Tourniquet
- \* Apply pressure to the wound
- \* Pack the wound to stop bleeding
- \* Placing a Israeli type bandage





Israeli Type Bandage



## How do we package our kits?

Each kit is vacuum sealed to ensure the product is sterile, waterproof, tamper resistant and convenient to open





guardiancprllc.com 920-229-4047